

Innovative Technology
For Your Medical Aesthetic Needs


COOLSHAPING

Extreme Cooling Lipolysis device



Treatment by cooling lipolysis technology

Coolshaping2 is a cooling lipolysis device. Combination of Cryo and Vacuum technology creates a safe and secure body-sliming effect. Compared to the existing device, the treatment area has been enlarged and the treatment time has been greatly reduced. Various types of 4 handpieces allow the treatment of local areas such as face.



Coolshaping's Revolutionary Change!

COOLSHAPING²

- ✓ Now, 360° All Area Cooling!
- ✓ Simultaneous operation of 4 handpieces!
- ✓ 4 Types of Handpiece for various treatment areas!
- ✓ Enlarged treatment area!
- ✓ Greatly reduced treatment time!

Treatment Mechanism



Start Cryo-cooling

Apply coolshaping applicator at the target area, start cryo-cooling to reach target temperature.



Immediate Cryo-reaction

The fat cells in the target area begin responding to cryo-crystallization.



30 Days Metabolism process

Over the following weeks and months, fat cells are disrupted and eliminated gradually by metabolism.



90 Days Eliminated cells

The fat cells are eliminated.

Advantages



Stable Cryo energy



Various applicators




Shot Downtime



Simultaneous operation



User convenient interface




Mini handpiece
*70 x 139 x 200 (mm)

Small handpiece
*89 x 139 x 119 (mm)

Medium handpiece
*93 x 188.7 x 96 (mm)

Large handpiece
*115 x 225.5 x 120 (mm)

Treatment Area



Before & After



Accessories

