



PULLING AND HOLDING THAT BEAUTY





BARBED CUTTING

Skin Quality

Lifting

croquis[®] BARBED CUTTING is used for tissue fixation and repositioning. With absorbable component and neocollagenesis, BARBED CUTTING type provides a younger and beautiful look.

Material : PDO (Polydioxanone)
Contains : 30ea / box (6ea / Al pouch * 5)
Thread Type : 

Product	Needle / Cannula		Thread		Barb
	Gauge	Length (mm)	USP	Length (mm)	
3D Type Barb (Needle)	21G	60	3 - 0	110	 (3D)
		90		150	
3D Type Barb (L-Cannula)	18G	100	2	150	
	19G	60	0	105	
		90		140	
	20G	60	2 - 0	105	
		90		140	