






RETURNING TO THAT YOUTH



MONO

Tightening ●●● | Volume ●○○○

croquis[®] MONO stimulates neocollagenesis, which not only improves your overall skin quality but also fully absorbs into the body over time.

Material : PDO (Polydioxanone)
Contains : 60ea / box (10ea / Al pouch * 6)
Thread Type : _____

Needle		Thread	
Gauge	Length (mm)	USP	Length (mm)
25G	60	5 - 0	90
	90		150
27G	60	6 - 0	90
	40		50
29G	50		70
	25		30
30G	40		50