

Innovative Technology
For Your Medical Aesthetic Needs

Hydra Beauty 2

Advanced Skin Care System



Advantage





Hydra Beauty 2

Advanced Skin Care System

Hydra Beauty 2 is an "Advanced Skin Care System" that protects the beauty of natural and healthy skin with customized treatments (exfoliation, pore, acne removal/brightening/moisturizing) for various skin concerns even for sensitive skin. It effectively removes skin waste, exfoliates, and supplies nutrients using a compact handpiece that is exclusively designed for Hydra Beauty 2.



Micro-Vibration

Bubble detox deep cleansing - Get rid of impurities and skin waste from pores with fine carbon bubbles created by a specialized tip / cream, Hypo Tip / Massager Tip - The combination of Hypo / Massager tip and micro-vibration increases solution penetration and induces the discharge of skin waste for clear and healthy skin.




Hydra Charcoal Detox



Hypo Tip



Massage Tip



Aqua Peeling

Hydra Beauty 2's exclusive solution with radial tip effectively dissolves dead skin waste, sebum, impurities and blackheads.



Small Tip



Large Tip



Body Tip



Diamond Peeling

Exfoliates dead skin waste and impurities by vacuum pressure and maximizes the solution absorption resulting in smoother skin.




ES-10



ES-20




ES-30




Electroporation

Electroporation technology uses high voltage electric pulses to temporarily open holes in the skin and effectively allows penetration of active ingredients. Safe and effective by penetrating various solutions into the dermis without damaging the skin cells.



Hydrogen Therapy

By applying high concentrated hydrogen into the skin cells converting the active oxygen in cells to water, preventing skin aging.



Solutions



S1 AHA (Alpha Hydroxy Acid)

Applies lactic Acid, a soluble solution (AHA), to dissolve bacteria, sebum and cosmetic residue within the pores.



S2 BHA (Beta Hydroxy Acid)

Applies salicylic acid, an oil controlling solution (BHA), to cleans the pore from sebum, black head in the T-zone.



S3 Hydration (Moisturizing component)

Permeates portulaca extract to provide vitamin B supplement, anti-oxidant effect, moisturization, melanin-control, and supply nourishing ingredients.

Electroporation Technology Calming Mask



Calming Mask, using Electroporation technology, opens a path between cells and allows penetration of active ingredients to dermal layer.



Hydra Beauty 2

Advanced Skin Care System



Head office: 120, Gieopdosi-ro, Jijeong-myeon, Wonju-si, Gangwon-do, Korea.
Tel: +82-33-760-1700(Rep.) Fax: +82-33-746-8716

Seoul office: #1701~1702, 145, Gasan digital 1-ro, Gaumcheon-gu, Seoul, Korea.
Tel: +82-2-514-9713(Rep.) Fax: +82-2-514-9716

E-mail: es@esglobal.co.kr / www.esglobal.co.kr