

TEOSYAL RHA® 2



For dynamic
smoothing of
moderate lines,
particularly glabella
and forehead
wrinkles, nasolabial
folds and for lip
reshaping¹

TEOSYAL RHA® 2

Consistent improvements over time

with only 30% of patients back to
pre-treatment score at 1 year¹

1. Rzany B, Converset-Vielthel S, Hartmann M, Lamouy JC, Ribé N, Sito G, Noize-Pin C. Efficacy and Safety of 3 New Resilient Hyaluronic Acid Fillers, Crosslinked With Decreased BQDE, for the Treatment of Dynamic Wrinkles: Results of an 18-Month, Randomized Controlled Trial Versus Anonymous Anonymous Comparators. Dermatol Surg. 2019 Oct;45(10):1304-1314.

TEOSYAL RHA® 2

TECHNICAL SPECIFICITIES^{1,2}



Hyaluronic Acid (HA) Concentration	23 mg/mL
Formula composition	Crosslinked HA
Degree of Modification (MoD) BDDE crosslinking	3.1%
Rheological properties	High Stretch Moderate Strength
Indications¹	Moderate wrinkles & lips: • Upper third • e.g. forehead • e.g. glabella
Estimated duration*	12-18 months ²
Lidocaine	0.3%
Packaging	2x 1 mL
Needle	30G x 1/2"

MoD: Degree of Modification. **HA:** Hyaluronic Acid.
BDDE: 1,4-Butanediol Diglycidyl Ether.

*Average duration depends on several factors: patients' skin type, wrinkle severity, injection type and volume

1. TEOSYAL RHA® 2 – Instructions For Use

2. Razany B, Conversat-Vielhe S, Hartmann M, Larmou JC, Ribé N, Sita G, Noze-Pin C. Efficacy and Safety of 3 New Resilient Hyaluronic Acid Fillers, Crosslinked With Decreased BDDE, for the Treatment of Dynamic Wrinkles: Results of an 18-Month, Randomized Controlled Trial Versus Already Available Comparators. Dermatol Surg. 2019 Oct;45(10):1304-1314



TEOSYAL RHA® 2

RHEOLOGICAL PARAMETERS^{1,2}




- ● ○ STRETCH
- ● ○ STRENGTH
- ● ● STRUCTURAL INTEGRITY

HA concentration: 23 mg/mL

RHA Collection

BY TEOXANE



1. Faivre J, Giller M, Tremblay E, Trévisio P, Bourdon F. Advanced Concepts in Rheology for the Evaluation of Hyaluronic Acid Based Soft Tissue Fillers. Dermatol Surg. 2021;47(5):e159-e167.
2. TEOSYAL RHA® 2 – Instructions For Use.

RHA® 1 = Teosyal® RHA 1, RHA® 2 = Teosyal® RHA 2, RHA® 3 = Teosyal® RHA 3, RHA® 4 = Teosyal® RHA 4



TEOSYAL RHA® 2

MDR INDICATIONS

Treatment options*



- ● ● ○ STRETCH
- ● ○ ○ STRENGTH
- ● ● ● STRUCTURAL INTEGRITY

HA concentration: 23 mg/mL



SUPERFICIAL LAYERS

FOREHEAD LINES	✓
FOREHEAD DEPRESSION*	✓
EYEBROWS	✓
GLABELLAR LINES	✓
TEMPLES	✓
CROWS' FEET	✓
NASOLABIAL FOLD	✓
MARIONETTE LINES	✓

