

TEOSYAL RHA® 3



For dynamic
filling of moderate
to severe folds,
particularly
nasolabial folds
and marionette
lines ¹

1. TEOSYAL RHA® 3 – Instructions For Use.



TEOSYAL RHA® 3


Significant improvements

In nasolabial folds at 6 months
($p > 0.001$)¹

1. Rzany B, Converset-Viethel S, Hartmann M, Lamouy JC, Ribé N, Sitt G, Noize-Pin C. Efficacy and Safety of 3 New Resilient Hyaluronic Acid Fillers, Crosslinked With Decreased BDDE, for the Treatment of Dynamic Wrinkles: Results of an 18-Month, Randomized Controlled Trial Versus Already Available Comparators. Dermatol Surg. 2019 Oct;45(10):1304-1314.

TEOSYAL RHA® 3

TECHNICAL SPECIFICITIES^{1,2}



*Average duration depends on several factors: patients' skin type, wrinkle severity, injection type and volume

1. TEOSYAL RHA® 3- Instructions For Use

2. Rzany B, Converset-Vielhail S, Hartmann M, Larrouy JC, Ribé N, Sita G, Noze-Pin C. Efficacy and Safety of 3 New Resilient Hyaluronic Acid Fillers, Crosslinked With Decreased BDDE, for the Treatment of Dynamic Wrinkles: Results of an 18-Month, Randomized Controlled Trial Versus Already Available Comparators. Dermatol Surg. 2019

TEOXANE'S RHA COLLECTION

Hyaluronic Acid (HA) Concentration	23 mg/mL
Formula composition	Crosslinked HA
Degree of Modification (MoD) BDDE crosslinking	3.6 %
Rheological properties	High Strength Good Stretch
Indications ¹	Severe wrinkles, philtrum and lips: e.g. Nasolabial folds e.g. Marionette lines
Estimated duration*	12 - 18 months ²
Lidocaine	0.3%
Packaging	2x 1 mL
Needle/cannula	27G x 1/2"

MoD: Degree of Modification. **HA:** Hyaluronic Acid.
BDDE: 1,4-Butanediol Diglycidyl Ether.



TEOSYAL RHA® 3

RHEOLOGICAL PARAMETERS^{1,2}




- ● ○ ○ STRETCH
- ● ● ○ STRENGTH
- ● ● ○ STRUCTURAL INTEGRITY

HA concentration: 23 mg/mL

RHA Collection

BY TEOXANE



1. Faivre J, Gallier M, Tremblais E, Trévisio P, Bourdon F. Advanced Concepts in Rheology for the Evaluation of Hyaluronic Acid-Based Soft Tissue Fillers. Dermatol Surg. 2021;47(5):e159-e167.
2. TEOSYAL RHA® 3 – Instructions For Use.

RHA® 1 = Teosyal® RHA 1, RHA® 2 = Teosyal® RHA 2, RHA® 3 = Teosyal® RHA 3, RHA® 4 = Teosyal® RHA 4



TEOSYAL RHA® 3


MDR INDICATIONS^{1,2}

Treatment options*



- ○ ○ STRETCH
- ● ○ STRENGTH
- ● ● STRUCTURAL INTEGRITY

HA concentration: 23 mg/mL



	 SUPERFICIAL LAYERS	 DEEP LAYERS
CHEEK	✓	
NASOLABIAL FOLD	✓	
PHILTRUM	✓	✓
ORAL COMMISSURES	✓	✓
LIPS	✓	✓
JAW LINE	✓	
MARIONETTE LINES	✓	
CHIN	✓	

*Average duration depends on several factors: patients' skin type, wrinkle severity, injection type and volume.

1. TEOXAL RHA® 3 - Instructions For Use

2. Rzany B, Converse-Verhel S, Hartmann M, Larrouy JC, Ribe N, Silo G, Noz-Pin C. Efficacy and Safety of 3 New Restylane Hyaluronic Acid Fillers, Crosslinked With Decreased BDDE, for the Treatment of Dynamic Wrinkles: Results of an 18-Month, Randomized Controlled Trial Versus Already Available Comparators. Dermatol Surg. 2019

