

TEOSYAL RHA KISS®



For dynamic
subtle lip
reshaping¹



TEOSYAL RHA® KISS

TECHNICAL SPECIFICITIES¹



- ● ● ○ STRETCH
- ● ○ ○ STRENGTH
- ● ● ● STRUCTURAL INTEGRITY

HA concentration: 23 mg/mL

Hyaluronic Acid (HA) Concentration	23 mg/mL
Formula composition	Crosslinked HMW HA
Degree of Modification (MoD) BDDE crosslinking	3.1 %
Rheological properties	High stretch Moderate strength
Indications ¹	Subtle lip reshaping, including philtrum and oral commissures
Estimated duration*	Up to 12 months
Lidocaine	0.3%
Packaging	2x 0.7 mL
Needle	30G x 1/2"

*Average duration depends on several factors: patients' skin type, wrinkle severity.

¹Injection type and volume

1. TEOSYAL RHA® KISS – Instructions For Use

2. Rzany B, Converset-Vielhél S, Hartmann M, Larouy JC, Ribe N, Sio G, Noze-Pin C. Efficacy and Safety of 3 New Resilient Hyaluronic Acid Fillers, Crosslinked With Decreased BDDE, for the Treatment of Dynamic Wrinkles: Results of an 18-Month, Randomized Controlled Trial Versus Already Available Comparators. Dermatol Surg. 2019

MoD: Degree of Modification. HA: Hyaluronic Acid.
BDDE: 1,4-Butanediol Diglycidyl Ether.



TEOSYAL RHA® KISS

MDR INDICATIONS¹

Treatment options*



- ● ● ○ STRETCH
- ● ● ○ STRENGTH
- ● ● ● STRUCTURAL INTEGRITY

HA concentration: 23 mg/mL



SUPERFICIAL LAYERS

PHILTRUM	✓
LIPS	✓
COMMISSURES	✓

